

5 Maintenance Tips for Your DELTA™ XRF Analyzer



Your Olympus handheld XRF analyzer is an investment, so here are some maintenance tips to help extend its life:

- 1. Inspect the Window:** The analyzer's window helps keep out dust and debris. Clean it daily and, if it breaks, [order a replacement](#).
- 2. Update the Alloy Database:** Get the most up-to-date alloy database for free by calling our technical support team at 877-225-8380 or [emailing us](#).
- 3. Run a Hardware Check:** Run the daily CalCheck procedure to confirm that the instrument hardware components (i.e., tube, detector, filters, and safety features) are operating properly.
- 4. Update the Batteries:** Order replacement batteries for your DELTA analyzer by [visiting our web store](#).
- 5. Upgrade the Software:** Keep your software up to date. [Upgrade it online](#) or call our technical support team at 877-225-8380.

Olympus Support Services

We have two U.S. service centers for repairs. When you're ready to replace your unit, our latest Vanta™ series analyzers offer improved speed, connectivity, and ruggedness. We offer trade-in discounts and [same-day shipping](#) on new instruments. To trade in your analyzer, call us at 800-225-8330 or send us an [email](#).